

BloodSTOP[®]


CONTROLS BLEEDING FAST



ADVANCED SURFACE WOUND CARE HEMOSTAT

- Accelerates clot formation, easily inserted into tight oral spaces
- Stops bleeding fast; upon contact with blood, turns into gel to seal the wound
- The gel prevents contamination and creates moist environment to promote healing

Size: 1.3cm x 5.0cm (24 pcs)



- Placed the dry BloodSTOP[®] iX on the area of interest (dry glove hand/instrument).
- Activate the strip with normal saline/blood. Request the patient, assistant and dentist not to touch the strip until it sets in about 3 minutes.
- Check the edge of the strip to ensure it has set and no longer sticky.



FACTUAL CONTENT AND IMAGERY SOURCED FROM **DR. GABRIEL SIA**

CONTACT US

☎ Holiday Ng +60 17-223 5520

✉ amertisbiomed@gmail.com

📍 ETHICS BIO SDN BHD
B-05-09 Kompleks Perindustrian EmHub,
Persiaran Surian, Seksyen 3, Taman Sains Selangor,
Kota Damansara ,47810, Petaling Jaya, Selangor.